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

ZEN 2.3 (blue edition)  
BIO Applications

# ZEN Blue主界面

The screenshot shows the ZEN Blue software interface with the following components labeled in Chinese:

- 标题栏** (Title Bar): Located at the top left, pointing to the window title "ZEN brain PH3 b-tub eng-GFP-EDF\_(c2+c3+c4+c5).czi - ZEN 2.3 lite".
- 菜单栏** (Menu Bar): Located at the top right, pointing to the menu items "File Edit View Acquisition Graphics Tools Window Help".
- 工具栏** (Toolbar): Located below the menu bar, pointing to the icons for Camera and Processing.
- 左侧工具栏** (Left Tool Bar): Located on the left side, pointing to the "Method" panel which lists various image processing functions like "Image Export", "Fuse Image Subset", and "Image Import".
- 右侧工具栏** (Right Tool Bar): Located on the right side, pointing to the "Images and Documents" panel which shows a list of open files.
- 图像显示区域** (Image Display Area): The central area pointing to the main image of a cell with a 100 μm scale bar.
- 图像控制区域** (Image Control Area): Located at the bottom, pointing to the "Display" panel which includes a histogram and various adjustment sliders for zoom, gamma, and white/black levels.
- 状态栏** (Status Bar): Located at the very bottom, pointing to the system information and performance metrics.

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






brain PH3 b-tub  
eng-GFP-EDF\_(c2+c3+c4+c5)\_c2.tif

brain PH3 b-tub  
eng-GFP-EDF\_(c2+c3+c4+c5)\_c3.tif

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The screenshot displays the ZEISS ZEN software interface with the following components and steps:

- Processing Panel (Left):** The 'Processing' button is highlighted with a red box and labeled '1'. Below it, the 'Batch' button is highlighted with a red box and labeled '2'. In the 'Batch Method' list, 'Image Export' is highlighted with a red box and labeled '3'.
- Batch Processing Panel (Middle):** The 'Add...' button is highlighted with a red box and labeled '4'. Below it, a table lists files to be processed:

S	Consistent	File Name	Size	Method	Output Name	Output Storage Path
5						

Select function first; Add Files by pressing the '+' button.

- Add Files Dialog (Bottom):** The 'Add files' dialog is open, showing a file list. The file '2017-0705-D1ER-bleach2-Create Image Subs...' is highlighted with a red box and labeled '5'. The 'Open (O)' button is highlighted with a red box and labeled '6'.

Method Parameters (Left):

- File type: JPEG File Interchange Format (JPEG)
- Quality: 100
- Resize: 100%
- Export to: C:\Users\ZCCHQIAO\Music

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The screenshot displays the 'Batch Processing' window in the Zeiss software. On the left, the 'Batch Method' list has 'Image Export' selected. Below it, the 'Method Parameters' panel is open, showing settings for 'Tagged Image File Format (TIFF)'. A red box highlights this panel, with an arrow pointing to it from the text '8、此处设置导出文件格式，和单个文件导出参数设置类似'. In the center, a table lists files for batch processing. A red arrow points to the first row of the table, with the text '7、选择其中一个文件设置导出文件参数'.

S	Consistent	File Name	Size	Method	Output Name	Output Storage Path
?		C:\Users\ZCCHQIAO\Downl...	6.36 MB	Image Export		C:\Users\ZCCHQIAO\Downloads\FRET
?		C:\Users\ZCCHQIAO\Downl...	1.67 MB	Image Export		C:\Users\ZCCHQIAO\Downloads\FRET
?		C:\Users\ZCCHQIAO\Downl...	29.15 MB	Image Export		C:\Users\ZCCHQIAO\Downloads\FRET
?		C:\Users\ZCCHQIAO\Downl...	1.99 MB	Image Export		C:\Users\ZCCHQIAO\Downloads\FRET
?		C:\Users\ZCCHQIAO\Downl...	1.78 MB	Image Export		C:\Users\ZCCHQIAO\Downloads\FRET
?		C:\Users\ZCCHQIAO\Downl...	6.42 MB	Image Export		C:\Users\ZCCHQIAO\Downloads\FRET
?		C:\Users\ZCCHQIAO\Downl...	1.85 MB	Image Export		C:\Users\ZCCHQIAO\Downloads\FRET

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